Phospho-Src (Tyr416) (C4) rabbit mAb PE conjugate

www.abwizbio.com

Support: info@abwizbio.com Order: sales@abwizbio.com

#2072

Store at: 2-8°C

For Research Use Only. Not For Use In Diagnostic Procedures.

Applications	Detection	Clonality	Isotype
Flow Cytometr	ry N/A	Monoclonal	Rabbit IgGk
Format:	PE		
Cross Reactivity:	Predicted to work with mouse, rat and other homologues.		
Formulation:	1X PBS, 0.09% NaN3, 0.2% BSA		
Preparation:	Protein A+G		
Reactivity:	Human,Mouse		
Recommended Usage:	For flow cytometric staining, the suggested µL of staining volume. It is recommended application.	d use of this reagent is 5 μL per m that the reagent be titrated for opt	iillion cells or 5 μL per 100 imal performance for each
Immunogen:	A synthetic phospho-peptide correspondin	g to residues surrounding Tyr416	of human phospho Src
Description:	phosphorylates specific tyrosine-protein kinase sic (c-sic, sic), is a non-receptor tyrosine kinase. Sic phosphorylates specific tyrosine residues in other proteins and is involved in signal transduction and is important in human tumor development and poor prognosis in cancers. There are nine enzymes (Fyn, Yes, Hck, Lck, Lyn, Blk, Frk, Yrk, and Fgr) closely-related to Src and are collectively called the Src family kinases (SFKs)(1). SFKs represent the largest family of nonreceptor tyrosine kinases and interact directly with receptor tyrosine kinases, G-protein-coupled receptors (GPCRs), steroid receptors, signals transducers, and activators of transcription and molecules involved in cell adhesion and migration. These events lead to proliferation, cell growth, differentiation, cell shape, motility, migration, angiogenesis, and survival. Phosphorylation at Tyr416 in the activation loop of the kinase domain upregulates phospho Src activity, while phosphorylation at Tyr527 in the carboxy-terminal tail decreases activity. Src has been strongly implicated in the development, maintenance, progression, and metastatic spread of several human cancers such as prostate, lung, breast, and colorectal. Structurally, SFKs are highly related to one another and contain conserved structural elements between family members. These elements include the N-terminal Src Homology 4 domain (SH4), the Src Homology 3 domain (SH3), Src Homology 2 domain (SH2), a linker sequence, the tyrosine kinase domain, and the C-terminal tail (2). The N-terminal domain, SH4, serves as a site for myristoylation and targets SFKs to the cytoplasmic membrane. The SH3 domain binds amino acid sequences rich in proline residues (3). This domain is critical for Src activity, intracellular localization, and the recruitment and binding of Src substrates. The SH2 domain binds to short motifs containing phospho-tyrosines. Together, the SH2 and SH3 domains cooperate in regulating SFKs catalytic activity		



References:

Events

(1) Src family kinases, key regulators of signal transduction. Parsons SJ, Parsons JT Oncogene. 2004 Oct 18; 23(48):7906-9.

(2) Src protein-tyrosine kinase structure and regulation. Roskoski R Jr Biochem Biophys Res Commun. 2004 Nov 26; 324(4):1155-64.

(3) Cohen GB, Ren R, Baltimore D. Modular binding domains in signal transduction proteins. Cell. 1995;80:237?248.



Flow cytometric analysis of U937 cells unstained untreated cells as negative control (blue) or stained untreated (red) or treated with IFN? IL4 and pervanadate (green) using phospho-Src (Tyr416) antibody SrcY416-C4 PE conjugate. Cat. #2072.

Phospho-Src (Y416) PE Coniudate

