

# SARS-CoV2 G10xA1 human neutralizing mAb

[www.abwizbio.com](http://www.abwizbio.com)

Support: [info@abwizbio.com](mailto:info@abwizbio.com)

Order: [sales@abwizbio.com](mailto:sales@abwizbio.com)

**Catalog:** #2571

**Store at:** -20°C

*For Research Use Only. Not For Use In Diagnostic Procedures.*

Applications	Detection	Clonality	Isotype
Functional Assay,ELISA	Anti-SARS-CoV-2 NP mAb	Monoclonal	Human IgG1k

**Format:** Unconjugated

**Cross Reactivity:** No

**Formulation:** 1X PBS

**Preparation:** Protein A

**Reactivity:** Other











**Recommended Usage:** A recombinant SARS-CoV-2 spike protein containing RBD.

**Immunogen:** A recombinant SARS-CoV-2 spike protein containing RBD.

**Description:**

**References:**

G10xA1 Surrogate Virus Neutralization Test IC50 (ng/mL)



N/A