Phospho-MCM2 (Ser139) (B12) rabbit mAb APC conjugate

www.abwizbio.com
Support: info@abwizbio.com
Order: sales@abwizbio.com

Catalog: #2434 **Store at:** 2-8ºC

For Research Use Only. Not For Use In Diagnostic Procedures.

Applications	Detection	Clonality	Isotype
Flow Cytometry	N/A	Monoclonal	Rabbit IgGk

Format: APC

Cross Reactivity: Predicted to work with mouse, rat and other homologues.

Formulation: 1X PBS, 0.09% NaN3, 0.2% BSA

Preparation: Protein A+G

Reactivity: Human,Mouse,Rat


Recommended Usage: For flow cytometric staining, the suggested use of this reagent is 5 µL per million cells or 5 µL per 100 µL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. See product image legends for additional information.

Immunogen: A synthetic phospho-peptide corresponding to residues surrounding Ser139 of human phospho MCM2

Description: The members of minichromosome maintenance (McM) protein family 2-7 were originally identified as a group of proteins essential for DNA replication (chromosomal maintenance (1,2). They share common sequence homology to each other in their nulceotide-binding domains and are distinct subgroup of the large AAA ATPase family, which are required for the initiation and elongation of DNA replication. It has been reported that Cdc7/Dbf4 phospohrylates MCM2 during G1/S cell cycle which coincides with the intiation of DNA replication (3,4)

References:

1. Bell S.P. and Dutta A. (2002) Annu Rev Biochem71:333-74.
2. Chong J.P. et al., (2995) Natur 375: 418-21.
3. Donovan S. et al., 1997) Proc Natl Acad Sci USA, 94:5611-6.
4. Tsuji T et al., (2006) Mol Biol Cell17:4459-72.



Flow cytometric analysis of C6 cells, untreated and unstained as negative control (blue) or untreated (red) or treated with staurosporine (green) and stained using Phospho-MCM2 (Ser139) antibody MCM2S139-B12 APC conjugate. Cat. #2434.

©2026 Abwiz Bio